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INTRODUCTION

This instrument is manufactured with the latest technology and requires no particular maintenance. PEQLAB certifies that this instrument was thoroughly inspected by the manufacturer, tested at the factory prior to shipment and found to meet all requirements defined by contract under which it is furnished. However, dimensions and other physical characteristics may vary.

The normal operating temperature should be between 0° and 40°C. Never use the instrument in a room with high humidity (>95 %) or at very low temperatures (condensation water!).

Connect the instrument only to an earthed power line. The required power source is indicated on the label at the back of the instrument. Do not cut and splice the power cord. When removing the power cord from the wall outlet, be sure to unplug by holding the plug attachment and not by pulling the cord. Do not hold the plug with wet hands.

DISTRIBUTOR

PEQLAB Biotechnologie GmbH
Carl-Thiersch-Str. 2b
91052 Erlangen
Germany

Tel (+49) (0) 9131 610 7020
Fax (+49) (0) 9131 610 7099
E-mail: info@peqlab.de

MANUFACTURER

Consort nv
Parklaan 36
B2300 Turnhout
Belgium

Tel (+32)(14)41 12 79
Fax (+32)(14)42 91 79
E-mail: info@Consort.be

WARRANTY

This instrument (excluding all accessories) is warranted against defective material and workmanship for a period of three years from the date of shipment. PEQLAB will repair all defective equipment returned during warranty without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. Please contact the PEQLAB technical service before returning a product for warranty repair! PEQLAB is not liable for consequential damages resulting from the use or handling of its products.

SERVICING

In the event of this instrument being returned for servicing, the owner is requested to remove the power supply lead and NOT to send the following items unless they are suspect:

- Manual
- Cables
- Accessories

If serious malfunctioning occurs, stop using the unit immediately and consult the PEQLAB technical service.

KEYBOARD

↑↓ = Button for entering a value or for selecting a function.

SET/ENTER = Button for programming the desired parameters.

RUN/STOP = Button to start or stop an experiment.

VERY IMPORTANT

This electrophoresis power supply is a high technology instrument available in several versions. As it is capable of giving dangerous voltage levels by which high power is involved, we suggest that you take a few moments to read this manual thoroughly. Although this instrument is equipped with all necessary safety features against abuse and other accidental failures, caution should be exercised when working with high voltage equipment.

- **Avoid to touch the outlets with any conducting object** and make sure there is a second person present for your safety in case of any severe electric shock.
- Never touch any part of the assembly (power supply, leads or tank) before having switched off.
- **Never manipulate with wet hands.**
- **Do not ground any of the outputs or the buffer in the tank.**
- Connect the outlets only to an insulated electrophoresis tank with safety cover.
- Never make any other connections, such as putting several power supplies in series or in parallel. In order to prevent electric shock, never open the back plate or remove the cover.
- Do not expose the unit to rain or any other liquid.
- **Do not spill liquid or insert metal objects into the unit.**
- Take care so that the power supply is not dropped to avoid damaging the cabinet which may defeat safeguards.
- If the unit has been dropped or the cabinet has been damaged, unplug it and have it checked by an authorised service technician to restore the safeguards.
- The fact that the unit operates satisfactorily does not imply that the unit is properly earthed or that it is completely safe. If there is any doubt about the effectiveness of the earthing of the unit, contact a qualified electrician.
- **Never block the ventilation holes** or place the unit in any enclosure unless proper ventilation is provided.
- Never place the unit near or over a radiator, heat register or stove.
- Avoid locations where the instrument is exposed directly to the sunlight.

SET-UP

1. Select SET-UP and press SET.
2. Select LANGUAGE and press SET.
3. Select the desired language and press SET.
4. Select CONTRAST and press SET.
5. Adjust the contrast of the display and press SET.
6. Select ALARM and press SET.
7. Decide whether or not the output voltage should shut down at a too low output current and press SET.
8. Select POWER FAIL and press SET.
9. Decide whether or not the instrument should proceed with the experiment for the remaining time after a power failure and press SET.
10. Select STAND-BY and press SET to return.

Select INFO and press SET to see general information about the electrophoresis power supply.

---- **STAND-BY** ----
SET-UP

---- **SET-UP** ----
LANGUAGE

CHOOSE LANGUAGE:
ENGLISH

---- **SET-UP** ----
CONTRAST

---- **SET-UP** ----
ALARM

CURRENT < 1 mA
ALARM: YES

---- **SET-UP** ----
POWER FAIL

POWER FAIL
DETECTION: YES

---- **SET-UP** ----
STAND-BY

---- **SET-UP** ----
INFO

METHOD PROGRAMMING

Method programming permits to store up to 9 different methods, each with 9 steps, of frequently used parameters in the non-volatile memory for future recall. Each step is able to recall a next one, providing a flexible multiple step function for special techniques. The method mode also permits to program a linear voltage gradient for any step provided the limiting current or power is not attained. Parameters can be changed temporarily without interrupting the run.

1. Verify that the instrument is switched off.
2. Connect the electrophoresis tank to the instrument.
3. Switch the instrument on.

If the timer is still running after switching on, the power supply will automatically proceed the experiment for the remaining time. Press STOP (within 10 seconds) to enter the stand-by mode to eventually change parameters. Press RUN to proceed with the experiment or press SET for 2 seconds to re-enter the stand-by mode.

4. Select the METHOD mode and press SET.
5. Select the desired method and press SET.
6. Select the desired step and press SET.
7. Decide whether or not to run a voltage gradient and press SET.
8. Enter the desired voltage limit and press SET (= the voltage end point in case of a gradient).
9. Enter the desired current limit and press SET (this limit is automatically set to maximum and cannot be changed in case of a voltage gradient).
10. Enter the desired power limit and press SET (this limit is automatically set to maximum and cannot be changed in case of a voltage gradient).
11. Enter the desired time-units (h or kVh) and press SET.
12. Enter the desired time and press SET.
13. Decide if the experiment should automatically proceed with a next step (resuming from point 6) or stop after this step. Enter the desired procedure and press SET.
14. Press RUN to start the experiment.

- At any time, you can return in the stand-by mode by holding the SET button for 2 seconds.
- At any time, you can avoid programming unnecessary parameters. To do so, simply press RUN instead of SET after entering the last relevant

**TIMER ACTIVE!
STOP NOW? 10s**

**---- STAND-BY ----
METHOD**

**ENTER 1|1
METHOD: 1**

**ENTER 1|1
STEP: 1**

**ENTER 1|1
GRADIENT: NO**

**ENTER 1|1
VOLTAGE: 400V**

**ENTER 1|1
CURRENT: 300mA**

**ENTER 1|1
POWER: 50W**

**ENTER 1|1
TIME UNIT: h**


**ENTER 1|1
TIME: 99:59h**

value and the instrument will immediately start the experiment keeping the other programmed parameters in memory.

- To ensure that the desired parameter is kept constant, we suggest to pre-set all other parameters to maximum.
- A blinking unit of the controlling parameter appears when the regulation is complete.
- The timer will always count down

ENTER	1 1
END METHOD:	YES

RUN	1 1	99:59h
400V	300mA	50W

MANUAL PROGRAMMING

Manual programming permits to set voltage, current, power and time limits for a simple routine electrophoresis run. Parameters can be changed temporarily without interrupting the run.

1. Verify that the instrument is switched off.
2. Connect the electrophoresis tank to the instrument.
3. Switch the instrument on.

If the timer is still running after switching on, the power supply will automatically proceed the experiment for the remaining time. Press STOP (within 10 seconds) to enter the stand-by mode to eventually change parameters. Press RUN to proceed with the experiment or press SET for 2 seconds to re-enter the stand-by mode.

4. Select the MANUAL mode and press SET.
5. Select the desired voltage limit and press SET.
6. Select the desired current limit and press SET.
7. Select the desired power limit and press SET.
8. Enter the desired time-units (hours or kilovolt hours) and press SET

TIMER ACTIVE!	10s
STOP NOW?	

----	STAND-BY	----
	MANUAL	

ENTER	
VOLTAGE:	400V

ENTER	
CURRENT:	300mA

ENTER	
POWER:	50W

9. Enter the desired time and press SET (pre-set to [--:--] if no timing is required)
10. Press RUN to start the experiment.

- *At any time, you can return in the stand-by mode by holding the SET button for 2 seconds.*
- *At any time, you can avoid programming unnecessary parameters. To do so, simply press RUN instead of SET after entering the last relevant value and the instrument will immediately start the experiment keeping the other programmed parameters in memory.*
- *To ensure that the desired parameter is kept constant, we suggest to pre-set all other parameters to maximum.*
- *A blinking unit of the controlling parameter appears when the regulation is complete.*
- *The timer will always count down. However, when no time limit has been programmed it will count up.*

ENTER TIME UNIT:	h
-----------------------------	----------

ENTER TIME:	99:59h
------------------------	---------------

RUN	1 1	99:59h
400V	300mA	50W

CHANGING PARAMETERS DURING A RUN

Only the parameters of the running step can be changed.

1. Press SET to change the parameters during the run without interrupting it.
2. Change the actual voltage limit and press SET.
3. Change the actual current limit and press SET.
4. Change the actual power limit, and press SET.
5. Change the actual time and press SET.
6. The power supply will now continue the run.

- *A gradient step can never be changed during a run.*
- *Temporarily changed parameters are not kept in memory. other programmed parameters in memory.*
- *The instrument will automatically return to the run mode when no button has been pressed during 10 seconds.*

SET	1 1	99:59h
400	300mA	50W

SET	1 1	99:59h
400V	300	50W

SET	1 1	99:59h
400V	300mA	50

SET	1 1	99:59h
400V	300mA	50W

DATA LOGGING

1. Select SET-UP and press SET.
2. Select DATA and press SET.
3. Select DURATION and press SET.
4. Select the desired time for the data-logging and press SET.
5. Select LOG and press SET.
6. Decide whether or not all data (voltage, current, power, time) during the next experiments should be logged and press SET.
7. Select SET-UP and press SET to return.

- *Starting a new experiment will erase the previous data in memory.*
- *Stored data can be transferred to a computer via the RS232 output using a software (www.consort.be).*
- *At any time, you can return in the stand-by mode by holding the SET button for 2 seconds.*

----	STAND-BY SET-UP	----
------	----------------------------	------

----	SET-UP DATA	----
------	------------------------	------

-----	DATA DURATION	-----
-------	--------------------------	-------

LOG DATA FOR:	10h
--------------------------	------------

-----	DATA LOG	-----
-------	---------------------	-------

LOG DATA:	YES
----------------------	------------

-----	DATA SET-UP	-----
-------	------------------------	-------

ALARMS

GROUND LEAKAGE

A dangerous ground leakage is present (check thoroughly the entire assembly). Your power supply is NOT defective! See Frequently Asked Questions at the end of this manual.

*** **ERROR** ***
GROUND LEAKAGE

SHORT-CIRCUIT

Connected tank is short-circuited or has a very low resistance (check thoroughly the complete assembly). Your power supply is NOT defective! See Frequently Asked Questions at the end of this manual.

*** **ERROR** ***
SHORT CIRCUIT

OVERLOAD

Connected tank has a too low resistance (check specifications). Your power supply is NOT defective! See Frequently Asked Questions at the end of this manual.

*** **ERROR** ***
OVERLOAD

CURRENT < 1mA

The output current is lower than 1 mA (switch the alarm off in the SET-UP menu). Your power supply is NOT defective! See Frequently Asked Questions at the end of this manual.

*** **ERROR** ***
CURRENT < 1mA

CALL SERVICE

Failure in the electronics. Contact the PEQLAB technical service (service@peqlab.de)

*** **ERROR** ***
CALL SERVICE

You may also try a general RESET: switch the instrument on while keeping the SET button pressed. Attention! All stored programs and data will be erased.

FREQUENTLY ASKED QUESTIONS

What are the relations between Voltage, Current, Power and Resistance?

Power (W) = Voltage (V) x Current (A)

Resistance (Ω) = Voltage (V) / Current (A)

How does a power supply react after pressing RUN?

The internal generator will start building up the high voltage at the output terminals while voltage and current are constantly measured and power calculated. When one of the pre-set parameters is exceeded, the generator stops and will keep that parameter constant.

How important is the resistance of an electrophoresis unit?

The resistance of an electrophoresis unit depends on its size, gel thickness, amount of buffer, buffer conductivity and temperature. This resistance will normally decrease in time due to a slowly increasing temperature. Electrophoresis units which have a resistance below the minimum load resistance of a power supply will trigger an alarm! Read the output voltage and current during a run to measure the resistance and use above formula to calculate the value.

How to keep a constant voltage during a run?

Program the desired voltage and a higher current and power than the maximum expected values:

Current > Voltage / Resistance

Power > Voltage x Current

How to keep a constant current during a run?

Program the desired current and a higher voltage and power than the maximum expected values:

Voltage > Current x Resistance

Power > Voltage x Current

How to keep a constant power during a run?

Program the desired power and a higher voltage and current than the maximum expected values:

Voltage > Current x Resistance

Current > Voltage / Resistance

Why are my output values different from those of a similar experiment?

Either your programmed parameters are not equal to those described or the resistance of your electrophoresis unit is different (see above). It cannot be due to e.g. another model of power supply as the relations between voltage, current, power and resistance are monitored in the same way by any instrument (the electrophysical laws are unchangeable and apply independent of the instrument used!).

What about connecting more than one unit to the same power supply?

The outlets being in parallel each electrophoresis unit will be supplied with exactly the same voltage. However, current and power may differ due to differences between them even when exactly the same model, gel, buffers, etc... are used. Therefore, it is recommended to run several electrophoresis units only in the constant voltage mode on the same power supply.

TECHNICAL DATA

Specifications	EV222	EV243	EV231
VOLTAGE	0...200 V	0...400 V	0...300 V
CURRENT	0...200 mA	0...300 mA	0...1000 mA
POWER	0...20 W	0...50 W	0...150 W
PARAMETER RANGE	1...100% of full scale	1...100% of full scale	1...100% of full scale
TIMER	0...99:59 h	0...99:59 h	0...99:59 h
VOLT-HOURS	0...99.99 kWh	0...99.99 kWh	0...99.99 kWh
DISPLAY	LCD, 2x16 characters	LCD, 2x16 characters	LCD, 2x16 characters
RESOLUTION	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W
PROGRAMS	-	9x9 set of parameters	9x9 set of parameters
OUTPUTS	3 in parallel, 4 mm sockets	3 in parallel, 4 mm sockets	4 in parallel, 4 mm sockets
MIN. LOAD RESISTANCE	15 Ω	30 Ω	10 Ω
NO LOAD DETECTION	on/off, programmable	on/off, programmable	on/off, programmable
GROUND LEAKAGE DETECTION	✓	✓	✓
OVERLOAD DETECTION	✓	✓	✓
COMPUTER CONTROL	✓	✓	✓
DATA-LOGGING	3600 values	3600 values	3600 values
RS232	9600 b/s	9600 b/s	9600 b/s
AMBIENT TEMPERATURE	0...40°C	0...40°C	0...40°C
RELATIVE HUMIDITY	0...95%, non condensing	0...95%, non condensing	0...95%, non condensing
POWER REQUIREMENTS	210...250 VAC, 50/60 Hz, 75 W	210...250 VAC, 50/60 Hz, 75 W	210...250 VAC, 50/60 Hz, 200 W
DIMENSIONS (WxDxH)	24x20x13 cm	24x20x13 cm	31x26x15 cm
WEIGHT	3 kg	3 kg	6 kg

Specifications	EV265	EV202	EV261
VOLTAGE	0...600 V	0...300 V	0...600 V
CURRENT	0...500 mA	0...2000 mA	0...1000 mA
POWER	0...150 W	0...300 W	0...300 W
PARAMETER RANGE	1...100% of full scale	1...100% of full scale	1...100% of full scale
TIMER	0...99:59 h	0...99:59 h	0...99:59 h
VOLT-HOURS	0...99.99 kWh	0...99.99 kWh	0...99.99 kWh
DISPLAY	LCD, 2x16 characters	LCD, 2x16 characters	LCD, 2x16 characters
RESOLUTION	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W
PROGRAMS	9x9 set of parameters	9x9 set of parameters	9x9 set of parameters
OUTPUTS	4 in parallel, 4 mm sockets	4 in parallel, 4 mm sockets	4 in parallel, 4 mm sockets
MIN. LOAD RESISTANCE	30 Ω	5 Ω	15 Ω
NO LOAD DETECTION	on/off, programmable	on/off, programmable	on/off, programmable
GROUND LEAKAGE DETECTION	✓	✓	✓
OVERLOAD DETECTION	✓	✓	✓
COMPUTER CONTROL	✓	✓	✓
DATA-LOGGING	3600 values	3600 values	3600 values
RS232	9600 b/s	9600 b/s	9600 b/s
AMBIENT TEMPERATURE	0...40°C	0...40°C	0...40°C
RELATIVE HUMIDITY	0...95%, non condensing	0...95%, non condensing	0...95%, non condensing
POWER REQUIREMENTS	210...250 VAC, 50/60 Hz, 200 W	210...250 VAC, 50/60 Hz, 360 W	210...250 VAC, 50/60 Hz, 360 W
DIMENSIONS (WxDxH)	31x26x15 cm	31x26x15 cm	31x26x15 cm
WEIGHT	6 kg	10 kg	10 kg

Specifications	EV215	EV232	EV233
VOLTAGE	0...1200 V	0...3000 V	0...3000 V
CURRENT	0...500 mA	0...150 mA	0...300 mA
POWER	0...300 W	0...150 W	0...300 W
PARAMETER RANGE	1...100% of full scale	1...100% of full scale	1...100% of full scale
TIMER	0...99:59 h	0...99:59 h	0...99:59 h
VOLT-HOURS	0...99.99 kWh	0...99.99 kWh	0...99.99 kWh
DISPLAY	LCD, 2x16 characters	LCD, 2x16 characters	LCD, 2x16 characters
RESOLUTION	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W	1 V, 1 mA, 1 W
PROGRAMS	9x9 set of parameters	9x9 set of parameters	9x9 set of parameters
OUTPUTS	4 in parallel, 4 mm sockets	4 in parallel, 4 mm sockets	4 in parallel, 4 mm sockets
MIN. LOAD RESISTANCE	70 Ω	600 Ω	300 Ω
NO LOAD DETECTION	on/off, programmable	on/off, programmable	on/off, programmable
GROUND LEAKAGE DETECTION	✓	✓	✓
OVERLOAD DETECTION	✓	✓	✓
COMPUTER CONTROL	✓	✓	✓
DATA-LOGGING	3600 values	3600 values	3600 values
RS232	9600 b/s	9600 b/s	9600 b/s
AMBIENT TEMPERATURE	0...40°C	0...40°C	0...40°C
RELATIVE HUMIDITY	0...95%, non condensing	0...95%, non condensing	0...95%, non condensing
POWER REQUIREMENTS	210...250 VAC, 50/60 Hz, 360 W	210...250 VAC, 50/60 Hz, 200 W	210...250 VAC, 50/60 Hz, 360 W
DIMENSIONS (WxDxH)	31x26x15 cm	31x26x15 cm	31x26x15 cm
WEIGHT	10 kg	6 kg	10 kg

Specifications	EV262
VOLTAGE	0...6000 V
CURRENT	0...150 mA
POWER	0...300 W
PARAMETER RANGE	1...100% of full scale
TIMER	0...99:59 h
VOLT-HOURS	0...99.99 kWh
DISPLAY	LCD, 2x16 characters
RESOLUTION	1 V, 1 mA, 1 W
PROGRAMS	9x9 set of parameters
OUTPUTS	4 in parallel, 4 mm sockets
MIN. LOAD RESISTANCE	1200 Ω
NO LOAD DETECTION	on/off, programmable
GROUND LEAKAGE DETECTION	✓
OVERLOAD DETECTION	✓
COMPUTER CONTROL	✓
DATA-LOGGING	3600 values
RS2_2	9600 b/s
AMBIENT TEMPERATURE	0...40°C
RELATIVE HUMIDITY	0...95%, non condensing
POWER REQUIREMENTS	210...250 VAC, 50/60 Hz, 360 W
DIMENSIONS (WxDxH)	31x26x15 cm
WEIGHT	10 kg

CERTIFICATES

PEQLAB WARRANTY CERTIFICATE

This instrument (excluding all accessories) is warranted against defective material and workman-ship for a period of three years from the date of shipment. PEQLAB will repair all defective equipment returned during warranty without charge, provided the equipment has been used under normal laboratory conditions and in accordance with the operating limitations and maintenance procedures in this instruction manual and when not having been subject to accident, alteration, misuse or abuse. Please contact the PEQLAB technical service before returning a product for war-ranty repair!

DECLARATION OF CONFORMITY

We declare under our sole responsibility that the product

Electrophoresis Power Supply content of the type numbers
EV202, EV215, EV231, EV232, EV233, EV243, EV261, EV262, EV26_

to which this declaration relates is in conformity with the following standards

EN61010

LOW VOLTAGE DIRECTIVE 73/23/EEG

EN50081-1

EN50082-1

EN60555-2

EMC DIRECTIVE 89/336/EEG

Turnhout, February 22, 2005

A handwritten signature in black ink, appearing to be 'F. Gans', is written over a horizontal line.

CERTIFICATE OF REGISTRATION

This is to certify that:

Consort NV

Turnhout

Has been assessed by QualityMasters in respect of its Quality Management System and found to comply with:

NBN EN ISO 9001:2000

Approval is hereby granted for registration providing the rules and conditions relating to certification are observed at all times.

Scope: Production of high-end equipment for electrochemistry and electrophoresis.

Original approval: 2005-08-29

Certificate expiry: 2010-08-29

Current certificate: 2005-09-17

Certificate number: BE-0211-317

Authorized signatory

For QualityMasters Belgium

This certificate remains the property of QualityMasters.

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PEQLAB Ltd.